




MODEL RCR4304A

CHROMOTOGRAPHY REFRIGERATOR



FRONT VIEW



SIDE VIEW

10/2005

REVCO
 275 Aiken Road
 Asheville, NC 28804
 TEL. 1-800-252-7100
 TEL. 828-658-2711 FAX 828-658-0363

MODEL RCR4304A

GENERAL SPECIFICATIONS

EXTERIOR - WHITE TEXTURED VINYL ON STEEL

WIDTH 1397mm (55") DEPTH 838mm (33") HEIGHT 1930mm (76")

INTERIOR - WHITE TEXTURED VINYL ON STEEL

WIDTH 1245mm (49") DEPTH 673mm (26.5") HEIGHT 1448mm (57")

INSULATION - FOAMED-IN-PLACE POLYURETHANE

GASKET - MAGNETIC

DOORS - (2) SLIDING GLASS, NO LOCKS

SHELVES - (8) EPOXY, ADJUSTABLE IN 51mm (2") INCREMENTS

CONTROL - (1) THERMOSTATIC, MAY BE SET FOR DESIRED TEMPERATURE

+2°C TO +10°C (36°F TO 50°F), FACTORY SET TO +4°C (39°F)

DEFROST - AUTOMATIC

COMPRESSOR - 1/3 H.P., R134a REFRIGERANT, NOMINAL BTU'S TO ROOM

AMBIENT - 690 WATTS (2355 BTU'S/HR)

ACCESS PORTS - (2) ONE ON EACH END

DECIBEL RATING - 73 dBA

CAPACITY - 1218 L (43 CU. FT.)

ELECTRICAL DATA - 120/60/1 19 AMPS

30 AMP PLUG, NEMA STYLE 5-30

(1) DUPLEX ELECTRICAL OUTLET LOCATED ON INSIDE BACK WALL

SHIPPING WEIGHT - 293 KG (645 LBS)

SHIPPING VOLUME - 3.4 CU. METERS (119 CU. FEET)