




MODEL RLR294A

LABORATORY REFRIGERATOR



FRONT VIEW



SIDE VIEW

3/2005

REVCO

REVCO
275 Aiken Road
Asheville, NC 28804
TEL. 1-800-252-7100
TEL. 828-658-2711 FAX 828-658-0363

MODEL RLR294A

GENERAL SPECIFICATIONS

EXTERIOR - WHITE TEXTURED VINYL ON STEEL

WIDTH 889mm (35") DEPTH 861mm (33.9") HEIGHT 1935mm (76.2")

INTERIOR - ALUMINUM, BAKED ENAMEL FINISH

WIDTH 813mm (32") DEPTH 660mm (26") HEIGHT 1654mm (65.1")

INSULATION - FOAMED-IN-PLACE POLYURETHANE

GASKET - MAGNETIC

DOOR - SOLID WITH KEY LOCK

SHELVES - (4) ZINC-PLATED, ADJUSTABLE IN 25mm (1") INCREMENTS

CASTERS - (4) WITH BRAKES

CONTROL - (1) THERMOSTATIC, MAY BE SET FOR DESIRED TEMPERATURE

+2°C TO +10°C (36°F TO 50°F), FACTORY SET TO +4°C (39°F)

DEFROST - AUTOMATIC

COMPRESSOR - 1/3 H.P., R134a REFRIGERANT, NOMINAL BTU'S TO ROOM

AMBIENT - 265 WATTS (903 BTU'S/HR)

DECIBEL RATING - 45 dBA

CAPACITY - 827 L (29.2 CU. FT.)

ELECTRICAL DATA - 120/60/1 5.7 AMPS

SHIPPING WEIGHT - 159 KG (350 LBS)

SHIPPING VOLUME - 2.4 CU. METERS (86 CU. FEET)