




# MODEL RLR4304A

## LABORATORY REFRIGERATOR



FRONT VIEW



SIDE VIEW

11/2003

# REVCO

REVCO  
275 Aiken Road  
Asheville, NC 28804  
TEL. 1-800-252-7100  
TEL. 828-658-2711 FAX 828-658-0363

# MODEL RLR4304A

## GENERAL SPECIFICATIONS

EXTERIOR - WHITE TEXTURED VINYL ON STEEL

WIDTH 1346mm (53") DEPTH 813mm (32") HEIGHT 1930mm (76")

INTERIOR - WHITE TEXTURED VINYL ON STEEL

WIDTH 1245mm (49") DEPTH 673mm (26.5") HEIGHT 1448mm (57")

INSULATION - FOAMED-IN-PLACE POLYURETHANE

GASKET - MAGNETIC

DOORS - (2) SLIDING GLASS, NO LOCKS

SHELVES - (8) EPOXY, ADJUSTABLE IN 51mm (2") INCREMENTS

CONTROL - THERMOSTATIC, MAY BE SET FOR DESIRED TEMPERATURE

+2°C TO +10°C (36°F TO 50°F), FACTORY SET TO +4°C (39°F)

DEFROST - AUTOMATIC

COMPRESSOR - 1/3 H.P., R134a REFRIGERANT, NOMINAL BTU'S TO ROOM

AMBIENT - 690 WATTS (2355 BTU'S/HR)

DECIBEL RATING - 73 dBA

CAPACITY - 1218 L (43 CU. FT.)

ELECTRICAL DATA - 120/60/1 9.2 AMPS

SHIPPING WEIGHT - 293 KG (645 LBS)

SHIPPING VOLUME - 3.4 CU. METERS (119 CU. FEET)